



OVOS1 mouse mAb

Catalog No	YP-mAb-08392
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	OVOS1
Protein Name	OVOS1
Immunogen	Synthesized peptide derived from human OVOS1 AA range: 61-111
Specificity	This antibody detects endogenous levels of OVOS1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:Is able to inhibit all four classes of proteinases by a unique 'trapping' mechanism.,similarity:Belongs to the protease inhibitor I39 (alpha-2-macroglobulin) family.,subunit:Homotetramer.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using OVOS1 mouse mAb

