




TLX2 mouse mAb

Catalog No	YP-mAb-08393
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TLX2 HOX11L1 NCX
Protein Name	TLX2
Immunogen	Synthesized peptide derived from human TLX2 AA range: 17-67
Specificity	This antibody detects endogenous levels of TLX2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Contains 1 homeobox DNA-binding domain.,
Background	This gene is a member of an orphan homeobox-containing transcription factor family. Studies of the mouse ortholog have shown that the encoded protein is crucial for the development of the enteric nervous system; in humans, loss-of-function may play a role in tumorigenesis of gastrointestinal stromal tumors. [provided by RefSeq, May 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TLX2 mouse mAb