




NFXL1 mouse mAb

Catalog No	YP-mAb-08399
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NFXL1 OZFP
Protein Name	NFXL1
Immunogen	Synthesized peptide derived from human NFXL1 AA range: 103-153
Specificity	This antibody detects endogenous levels of NFXL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the NFX1 family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 10 NF-X1-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NFXL1 mouse mAb