




UAP1L mouse mAb

Catalog No	YP-mAb-08404
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UAP1L1
Protein Name	UAP1L
Immunogen	Synthesized peptide derived from human UAP1L AA range: 287-337
Specificity	This antibody detects endogenous levels of UAP1L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	56kD
Cell Pathway	cytosol,
Tissue Specificity	
Function	similarity:Belongs to the UDPGP type 1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night