




U17L6 mouse mAb

Catalog No	YP-mAb-08429
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	USP17L6 USP17C USP17D USP17N
Protein Name	U17L6
Immunogen	Synthesized peptide derived from human U17L6 AA range: 207-257
Specificity	This antibody detects endogenous levels of U17L6 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.;caution:The RS447 megasatellite DNA is a highly polymorphic conserved tandem repetitive sequence which contains a copy of USP17. It is present with an interindividual variation in copy number and between 20 to 103 copies can be found in the genome both on chromosome 4 and chromosome 8.;function:Seems to be non-functional in the regulation of apoptosis. Cleaves ubiquitin fusion protein substrates.;similarity:Belongs to the peptidase C19 family. USP17 subfamily.;
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using U17L6 mouse mAb