




5NT3L mouse mAb

Catalog No	YP-mAb-08456
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	NT5C3L
Protein Name	5NT3L
Immunogen	Synthesized peptide derived from human 5NT3L AA range: 23-73
Specificity	This antibody detects endogenous levels of 5NT3L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	catalytic activity:A 5'-ribonucleotide + H(2)O = a ribonucleoside + phosphate.,function:Can act both as nucleotidase and as phosphotransferase.,similarity:Belongs to the pyrimidine 5'-nucleotidase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using 5NT3L mouse mAb