




MSMP mouse mAb

Catalog No	YP-mAb-08463
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MSMP PSMP
Protein Name	MSMP
Immunogen	Synthesized peptide derived from human MSMP AA range: 58-108
Specificity	This antibody detects endogenous levels of MSMP at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	Detected in prostate epithelium (at protein level) (PubMed:17338636). Detected in trachea and testis (PubMed:17338636). Highly expressed in benign prostatic hyperplasia and in some prostate cancers, and can also be detected in breast tumor tissue (PubMed:17338636, PubMed:24442440).
Function	similarity:Belongs to the beta-microseminoprotein family.,tissue specificity:Detected in prostate epithelium (at protein level). Detected in trachea and testis. Highly expressed in the human prostate cancer cell line PC-3.,
Background	This gene encodes a member of the beta-microseminoprotein family. Members of this protein family contain ten conserved cysteine residues that form intra-molecular disulfide bonds. The encoded protein may play a role in prostate cancer tumorigenesis. [provided by RefSeq, Jan 2011],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using MSMP mouse mAb