



UBR3 mouse mAb

Catalog No	YP-mAb-08500
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UBR3 KIAA2024 ZNF650
Protein Name	UBR3
Immunogen	Synthesized peptide derived from human UBR3 AA range: 1591-1641
Specificity	This antibody detects endogenous levels of UBR3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	function:E3 ubiquitin-protein ligase which is a component of the N-end rule pathway. Does not bind to proteins bearing specific N-terminal residues that are destabilizing according to the N-end rule, leading to their ubiquitination and subsequent degradation.,pathway:Protein modification; protein ubiquitination.,sequence caution:Translated as Gln.,similarity:Belongs to the UBR1 family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 1 UBR-type zinc finger.,subunit:Interacts with UBE2A and UBE2B.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images