



COLL1 mouse mAb

Catalog No	YP-mAb-08502
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CLPSL1 C6orf127
Protein Name	COLL1
Immunogen	Synthesized peptide derived from human COLL1 AA range: 35-85
Specificity	This antibody detects endogenous levels of COLL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	Exclusively expressed in epididymis, in the corpus region.
Function	similarity:Belongs to the colipase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using COLL1 mouse mAb

