




ORML1 mouse mAb

Catalog No	YP-mAb-08521
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ORMDL1 HSPC202
Protein Name	ORML1
Immunogen	Synthesized peptide derived from human ORML1 AA range: 22-72
Specificity	This antibody detects endogenous levels of ORML1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Expressed in adult and fetal heart, brain, lung, liver, skeletal muscle and kidney. Expressed in adult pancreas and placenta and in fetal spleen abd thymus. Expressed at intermediate level in pancreas, placenta and brain but low in skeletal muscle and lung.
Function	similarity:Belongs to the ORM family.,tissue specificity:Widely expressed. Expressed in adult and fetal heart, brain, lung, liver, skeletal muscle and kidney. Expressed in adult pancreas and placenta and in fetal spleen abd thymus. Expressed at intermediate level in pancreas, placenta and brain but low in skeletal muscle and lung.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ORML1 mouse mAb