




FRG2 mouse mAb

Catalog No	YP-mAb-08537
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FRG2 FRG2A
Protein Name	FRG2
Immunogen	Synthesized peptide derived from human FRG2 AA range: 104-154
Specificity	This antibody detects endogenous levels of FRG2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Expression is undetectable in all tissues tested except for differentiating myoblasts of FSHD patients, which display low, yet distinct levels of expression, partly from FRG2, but predominantly originating from its homolog on chromosome 10.
Function	similarity:Belongs to the FRG2 family.,tissue specificity:Expression is undetectable in all tissues tested except for differentiating myoblasts of FSHD patients, which display low, yet distinct levels of expression, partly from FRG2, but predominantly originating from its homolog on chromosome 10.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FRG2 mouse mAb