




TULP4 mouse mAb

Catalog No	YP-mAb-08542
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TULP4 KIAA1397 TUBL4 TUSP
Protein Name	TULP4
Immunogen	Synthesized peptide derived from human TULP4 AA range: 798-848
Specificity	This antibody detects endogenous levels of TULP4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm.
Tissue Specificity	Expressed mainly in the brain, skeletal muscle, testis and kidney.
Function	domain:The SOCS box domain mediates the interaction with the Elongin BC complex, an adapter module in different E3 ubiquitin ligase complexes.,function:May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the TUB family.,similarity:Contains 1 SOCS box domain.,similarity:Contains 3 WD repeats.,tissue specificity:Expressed mainly in the brain, skeletal muscle, testis and kidney.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TULP4 mouse mAb