




ANKY1 mouse mAb

Catalog No	YP-mAb-08641
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ANKMY1 TSAL1 ZMYND13
Protein Name	ANKY1
Immunogen	Synthesized peptide derived from human ANKY1 AA range: 886-936
Specificity	This antibody detects endogenous levels of ANKY1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 MYND-type zinc finger.,similarity:Contains 3 MORN repeats.,similarity:Contains 7 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ANKY1 mouse mAb