




LRC4B mouse mAb

Catalog No	YP-mAb-08644
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	LRRC4B LRIG4
Protein Name	LRC4B
Immunogen	Synthesized peptide derived from human LRC4B AA range: 22-72
Specificity	This antibody detects endogenous levels of LRC4B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane; Single-pass membrane protein. Cell junction, synapse, presynaptic cell membrane .
Tissue Specificity	
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 9 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRC4B mouse mAb