




CQ097 mouse mAb

Catalog No	YP-mAb-08665
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C17orf97
Protein Name	CQ097
Immunogen	Synthesized peptide derived from human CQ097 AA range: 328-378
Specificity	This antibody detects endogenous levels of CQ097 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CQ097 mouse mAb