




YIPF4 mouse mAb

Catalog No	YP-mAb-08704
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	YIPF4 Nbla11189
Protein Name	YIPF4
Immunogen	Synthesized peptide derived from human YIPF4 AA range: 59-109
Specificity	This antibody detects endogenous levels of YIPF4 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in keratinocytes (at protein level).
Function	similarity:Belongs to the YIP1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using YIPF4 mouse mAb