




DJC13 mouse mAb

Catalog No	YP-mAb-08798
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DNAJC13 KIAA0678 RME8
Protein Name	DJC13
Immunogen	Synthesized peptide derived from human DJC13 AA range: 1576-1626
Specificity	This antibody detects endogenous levels of DJC13 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Early endosome . Early endosome membrane ; Peripheral membrane protein . Endosome membrane .
Tissue Specificity	
Function	sequence caution:Contaminating sequence. Potential poly-A sequence starting in position 748.,sequence caution:Contaminating sequence. Potential poly-A sequence.,sequence caution:Translated as Glu.,similarity:Contains 1 J domain.,
Background	This gene encodes a member of the Dnaj protein family whose members act as co-chaperones of a partner heat-shock protein by binding to the latter and stimulating ATP hydrolysis. The encoded protein associates with the heat-shock protein Hsc70 and plays a role in clathrin-mediated endocytosis. It may also be involved in post-endocytic transport mechanisms via its associations with other proteins, including the sorting nexin SNX1. Mutations in this gene are associated with Parkinson's disease. [provided by RefSeq, Jun 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DJC13 mouse mAb