




P3IP1 mouse mAb

Catalog No	YP-mAb-08805
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	PIK3IP1 HGFL
Protein Name	P3IP1
Immunogen	Synthesized peptide derived from human P3IP1 AA range: 107-157
Specificity	This antibody detects endogenous levels of P3IP1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	similarity:Contains 1 kringle domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using P3IP1 mouse mAb