




TRAIP mouse mAb

Catalog No	YP-mAb-08822
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRAIP RNF206 TRIP
Protein Name	TRAIP
Immunogen	Synthesized peptide derived from human TRAIP AA range: 399-449
Specificity	This antibody detects endogenous levels of TRAIP at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleoplasm . Nucleus, nucleolus . Chromosome . Cytoplasm. Cytoplasm, perinuclear region . In the nucleus, found in close proximity to PCNA, suggesting localization at replication foci (PubMed:26595769). Localizes to DNA damage sites in response to replication stress (PubMed:26781088, PubMed:26595769, PubMed:26711499). .
Tissue Specificity	
Function	function:Inhibits activation of NF-kappa-B mediated by TNF.,similarity:Contains 1 RING-type zinc finger.,subunit:Binds TRAF1 and TRAF2 and is part of the receptor-TRAF signaling complex (By similarity). Interacts with CYLD.,
Background	This gene encodes a protein that contains an N-terminal RING finger motif and a putative coiled-coil domain. A similar murine protein interacts with TNFR-associated factor 1 (TRAF1), TNFR-associated factor 2 (TRAF2), and cylindromatosis. The interaction with TRAF2 inhibits TRAF2-mediated nuclear factor kappa-B, subunit 1 activation that is required for cell activation and protection against apoptosis. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TRAIP mouse mAb