




DTWD1 mouse mAb

Catalog No	YP-mAb-08828
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	DTWD1 MDS009
Protein Name	DTWD1
Immunogen	Synthesized peptide derived from human DTWD1 AA range: 31-81
Specificity	This antibody detects endogenous levels of DTWD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	domain:Contains 1 DXTW motif.,similarity:Belongs to the DTW family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DTWD1 mouse mAb