




ZN276 mouse mAb

Catalog No	YP-mAb-08846
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZNF276 ZFP276 ZNF477
Protein Name	ZN276
Immunogen	Synthesized peptide derived from human ZN276 AA range: 325-375
Specificity	This antibody detects endogenous levels of ZN276 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Chromosome, centromere, kinetochore .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Contains 1 ZAD (zinc-finger associated) domain.,similarity:Contains 5 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZN276 mouse mAb