




DNJA2 mouse mAb

Catalog No	YP-mAb-08871
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	DNAJA2 CPR3 HIRIP4
Protein Name	DNJA2
Immunogen	Synthesized peptide derived from human DNJA2 AA range: 309-359
Specificity	This antibody detects endogenous levels of DNJA2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Lipid-anchor .
Tissue Specificity	
Function	function:Co-chaperone of Hsc70.,similarity:Contains 1 CR-type zinc finger.,similarity:Contains 1 J domain.,
Background	The protein encoded by this gene belongs to the evolutionarily conserved DNAJ/HSP40 family of proteins, which regulate molecular chaperone activity by stimulating ATPase activity. DNAJ proteins may have up to 3 distinct domains: a conserved 70-amino acid J domain, usually at the N terminus; a glycine/phenylalanine (G/F)-rich region; and a cysteine-rich domain containing 4 motifs resembling a zinc finger domain. The product of this gene works as a cochaperone of Hsp70s in protein folding and mitochondrial protein import in vitro. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DNJA2 mouse mAb