




RGL3 mouse mAb

Catalog No	YP-mAb-08932
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RGL3
Protein Name	RGL3
Immunogen	Synthesized peptide derived from human RGL3 AA range: 601-651
Specificity	This antibody detects endogenous levels of RGL3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	intracellular,
Tissue Specificity	
Function	domain:The Ras-associating domain plays a central role in the activation of Ral-A GDP/GTP exchange activity.,function:Guanine nucleotide exchange factor (GEF) for Ral-A. Potential effector of GTPase HRas and Ras-related protein M-Ras. Negatively regulates Elk-1-dependent gene induction downstream of HRas and MEKK1.,similarity:Contains 1 N-terminal Ras-GEF domain.,similarity:Contains 1 Ras-associating domain.,similarity:Contains 1 Ras-GEF domain.,subunit:Interacts with GTP-bound forms of RIT1, HRAS and MRAS.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RGL3 mouse mAb