




FA73B mouse mAb

Catalog No	YP-mAb-08962
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM73B C9orf54 PSEC0112
Protein Name	FA73B
Immunogen	Synthesized peptide derived from human FA73B AA range: 6-56
Specificity	This antibody detects endogenous levels of FA73B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion outer membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the FAM73 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA73B mouse mAb