




GG6L9 mouse mAb

Catalog No	YP-mAb-08964
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GOLGA6L9
Protein Name	GG6L9
Immunogen	Synthesized peptide derived from human GG6L9 AA range: 151-201
Specificity	This antibody detects endogenous levels of GG6L9 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	atus,
Tissue Specificity	
Function	similarity:Belongs to the GOLGA6 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GG6L9 mouse mAb