




CQ099 mouse mAb

Catalog No	YP-mAb-09012
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C17orf99 UNQ464/PRO809
Protein Name	CQ099
Immunogen	Synthesized peptide derived from human CQ099 AA range: 90-140
Specificity	This antibody detects endogenous levels of CQ099 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in fetal liver and bone marrow (PubMed:28978694). Expressed in peripheral blood lymphocyte B cells (PubMed:28978694).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CQ099 mouse mAb