




CHUR mouse mAb

Catalog No	YP-mAb-09019
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CHURC1 C14orf52 CHCH My015
Protein Name	CHUR
Immunogen	Synthesized peptide derived from human CHUR AA range: 51-101
Specificity	This antibody detects endogenous levels of CHUR at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	function:Transcriptional activator that mediates FGF signaling during neural development. Plays a role in the regulation of cell movement (By similarity). Does not bind DNA by itself.,mass spectrometry:The mass difference may be explained by bound zinc ions PubMed:17610897,similarity:Belongs to the Churchill family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CHUR mouse mAb