




S35G1 mouse mAb

Catalog No	YP-mAb-09037
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SLC35G1 C10orf60 TMEM20
Protein Name	S35G1
Immunogen	Synthesized peptide derived from human S35G1 AA range: 14-64
Specificity	This antibody detects endogenous levels of S35G1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Translocates from the endoplasmic reticulum to the cell membrane in response to a depletion of intracellular calcium and is detected at punctae corresponding to junctions between the endoplasmic reticulum and the cell membrane.
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Contains 2 DUF6 domains.,
Background	This gene encodes a transmembrane protein which is a member of the drug/metabolite transporter protein superfamily. The encoded protein may play a role in the regulation of calcium levels inside the cell. [provided by RefSeq, Sep 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using S35G1 mouse mAb