




TMPS7 mouse mAb

Catalog No	YP-mAb-09068
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TMPRSS7
Protein Name	TMPS7
Immunogen	Synthesized peptide derived from human TMPS7 AA range: 87-137
Specificity	This antibody detects endogenous levels of TMPS7 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Single-pass type II membrane protein .
Tissue Specificity	Expressed in brain, ovary, testis, salivary gland, trachea and lung.
Function	function:Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position.,PTM:N-glycosylated.,sequence caution:Differs at the N-terminus for unknown reasons.,sequence caution:Translation N-terminally extended.,similarity:Belongs to the peptidase S1 family.,similarity:Contains 1 peptidase S1 domain.,similarity:Contains 2 CUB domains.,similarity:Contains 2 LDL-receptor class A domains.,tissue specificity:Expressed in brain, ovary, testis, salivary gland, trachea and lung.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TMPS7 mouse mAb