




## ZSC29 mouse mAb

|                           |                                                                                                                                                                                                                        |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                | YP-mAb-09092                                                                                                                                                                                                           |
| Isotype                   | IgG                                                                                                                                                                                                                    |
| Reactivity                | Human;Rat;Mouse;                                                                                                                                                                                                       |
| Applications              | WB                                                                                                                                                                                                                     |
| Gene Name                 | ZSCAN29 ZNF690                                                                                                                                                                                                         |
| Protein Name              | ZSC29                                                                                                                                                                                                                  |
| Immunogen                 | Synthesized peptide derived from human ZSC29 AA range: 409-459                                                                                                                                                         |
| Specificity               | This antibody detects endogenous levels of ZSC29 at Human                                                                                                                                                              |
| Formulation               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                |
| Source                    | Monoclonal, Mouse,IgG                                                                                                                                                                                                  |
| Purification              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                   |
| Dilution                  | WB 1:500-1:2000                                                                                                                                                                                                        |
| Concentration             | 1 mg/ml                                                                                                                                                                                                                |
| Purity                    | ≥90%                                                                                                                                                                                                                   |
| Storage Stability         | -20°C/1 year                                                                                                                                                                                                           |
| Synonyms                  |                                                                                                                                                                                                                        |
| Observed Band             |                                                                                                                                                                                                                        |
| Cell Pathway              | Nucleus .                                                                                                                                                                                                              |
| Tissue Specificity        |                                                                                                                                                                                                                        |
| Function                  | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 6 C2H2-type zinc fingers., |
| Background                |                                                                                                                                                                                                                        |
| matters needing attention | Avoid repeated freezing and thawing!                                                                                                                                                                                   |
| Usage suggestions         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                      |



## Products Images



Western Blot analysis of various cells using ZSC29 mouse mAb