




ACRC mouse mAb

Catalog No	YP-mAb-09098
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ACRC
Protein Name	ACRC
Immunogen	Synthesized peptide derived from human ACRC AA range: 444-494
Specificity	This antibody detects endogenous levels of ACRC at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Nucleus, PML body . Chromosome . Co-localizes with SUMO2 at PML bodies in all interphase cells (PubMed:30914427). Localizes on condensed chromosomes in spermatocytes in G2 and M during meiotic prophase (By similarity). .
Tissue Specificity	Expressed in germ cells of the testis (at protein level) (PubMed:27718356). Detected in skeletal muscle, liver, kidney, pancreas, heart, lung and brain (PubMed:11714101). Expressed throughout spermatogenesis, from spermatogonia to elongated spermatids, in normal adult testis (at protein level) (PubMed:33963445).
Function	similarity:Belongs to the serine-aspartate repeat-containing protein (SDr) family.,tissue specificity:Detected in skeletal muscle, liver, kidney, pancreas, heart, lung and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ACRC mouse mAb