




CN132 mouse mAb

Catalog No	YP-mAb-09152
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C14orf132
Protein Name	CN132
Immunogen	Synthesized peptide derived from human CN132 AA range: 58-108
Specificity	This antibody detects endogenous levels of CN132 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CN132 mouse mAb