



LUZP4 mouse mAb

Catalog No	YP-mAb-09366
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LUZP4
Protein Name	Leucine zipper protein 4 (Cancer/testis antigen 28) (CT-28) (CT28) (Tumor antigen HOM-TES-85)
Immunogen	Synthesized peptide derived from human LUZP4 AA range: 42-92
Specificity	This antibody detects endogenous levels of human LUZP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serumby affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	34kD
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm . In nuclear speckles. Relocalizes to the cytoplasm during cell division. .
Tissue Specificity	Expressed specifically in testis. Also expressed in a wide variety of cancer types, but particularly high levels of expression observed in melanoma cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

