



ASCL4 mouse mAb

Catalog No	YP-mAb-09367
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ASCL4 BHLHA44 HASH4
Protein Name	Achaete-scute homolog 4 (ASH-4) (hASH4) (Achaete-scute-like protein 4) (Class A basic helix-loop-helix protein 44) (bHLHa44)
Immunogen	Synthesized peptide derived from human ASCL4 AA range: 81-131
Specificity	This antibody detects endogenous levels of human ASCL4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	19kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in skin. 7-fold higher expression in fetal skin than in adult skin. Weak expression also detected in fetal lung, aorta and brain, and in adult stomach, kidney, ovary and breast.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

