



LPP60 mouse mAb

Catalog No	YP-mAb-09615
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ASPG C14orf76
Protein Name	60 kDa lysophospholipase (EC 3.1.1.5) [Includes: L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase); Platelet-activating factor acetylhydrolase (PAF acetylhydrolase) (EC 3.1.1.47)]
Immunogen	Synthesized peptide derived from human LPP60 AA range: 112-162
Specificity	This antibody detects endogenous levels of human LPP60
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	63kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

