



AMFR mouse mAb

Catalog No	YP-mAb-10091
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AMFR RNF45
Protein Name	E3 ubiquitin-protein ligase AMFR (EC 6.3.2.-) (Autocrine motility factor receptor, isoform 2) (AMF receptor, isoform 2) (RING finger protein 45) (gp78)
Immunogen	Synthesized peptide derived from human AMFR AA range: 572-622
Specificity	This antibody detects endogenous levels of human AMFR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	71KD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Palmitoylation promotes localization to the peripheral endoplasmic reticulum. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images