



FO XK2 mouse mAb

Catalog No	YP-mAb-10102
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FO XK2 ILF ILF1
Protein Name	Forkhead box protein K2 (Cellular transcription factor ILF-1) (FO XK1) (Interleukin enhancer-binding factor 1)
Immunogen	Synthesized peptide derived from human FO XK2 AA range: 331-381
Specificity	This antibody detects endogenous levels of human FO XK2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	73kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Expressed in both lymphoid and non-lymphoid cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images