



IVD mouse mAb

Catalog No	YP-mAb-10148
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IVD
Protein Name	Isovaleryl-CoA dehydrogenase, mitochondrial (IVD) (EC 1.3.8.4)
Immunogen	Synthesized peptide derived from human IVD AA range: 325-375
Specificity	This antibody detects endogenous levels of human IVD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	47kD
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

