



Tubulin alpha mouse mAb(Yeast)

Catalog No	YP-mAb-10304
Isotype	IgG
Reactivity	Yeast
Applications	WB
Gene Name	
Protein Name	Tubulin alpha
Immunogen	Purified recombinant Full length Tubulin alpha (Yeast)
Specificity	This antibody detects endogenous levels of Tubulin alpha (Yeast)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tubulin alpha
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using Tubulin alpha mouse mAb(Yeast)

