



CSH1 Monoclonal Antibody

Catalog No	YP-mAb-10622
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CSH1; CSH2
Protein Name	chorionic somatomammotropin hormone 1 (placental lactogen)
Immunogen	Synthetic peptide from human protein at AA range: 161-210
Specificity	The antibody detects endogenous CSH1 protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSH1; CSH2
Observed Band	25kD
Cell Pathway	extracellular region,endoplasmic reticulum,vesicle,membrane-bounded vesicle,
Tissue Specificity	Placenta,Uterus,
Function	function:Similar to that of somatotropin.,miscellaneous:The sequence of CSH1 is shown.,similarity:Belongs to the somatotropin/prolactin family.,
Background	The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones and plays an important role in growth control. The gene is located at the growth hormone locus on chromosome 17 along with four other related genes in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. Although the five genes share a remarkably high degree of sequence identity, they are expressed selectively in different tissues. Alternative splicing generates additional isoforms of each of the five growth hormones, leading to further diversity and potential for specialization. This particular family member is expressed mainly in the placenta and utilizes multiple transcription initiation sites. Expression of the identical mature proteins for chorionic somatomammotropin hormones 1 and 2 is upregulated during development, al


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CSH1 Monoclonal Antibody