




Tmem72 Monoclonal Antibody

Catalog No	YP-mAb-10824
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM72 C10orf127 KSP37
Protein Name	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Immunogen	Synthesized peptide derived from human Tmem72 Monoclonal
Specificity	This antibody detects endogenous levels of Tmem72.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Observed Band	30kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep 2016],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



1 mouse-liver

2 CAC02

Western Blot analysis of various cells using Tmem72
Monoclonal Antibody