




TCEAL7 Polyclona Antibody

Catalog No	YP-mAb-10892
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TCEAL7
Protein Name	TCEAL7
Immunogen	Synthesized peptide derived from human TCEAL7 AA range: 34-84
Specificity	This antibody detects endogenous levels of human TCEAL7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Transcription elongation factor A protein-like 7 (TCEA-like protein 7;Transcription elongation factor S-II protein-like 7)
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in normal and fetal brain tissues, and weakly expressed in uterus and ovary. Down-regulated in epithelial ovarian, cervical, prostate, breast, brain and lung cancer cell lines and in brain and ovarian tumors.
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TCEAL7
Polyclonal Antibody