




F198B mouse mAb

Catalog No	YP-mAb-10934
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FAM198B C4orf18 ENED AD021 UNQ2512/PRO6001
Protein Name	F198B
Immunogen	Synthesized peptide derived from human F198B AA range: 372-422
Specificity	This antibody detects endogenous levels of F198B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using F198B mouse mAb