




GLYR1 mouse mAb

Catalog No	YP-mAb-10980
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	GLYR1 HIBDL NP60
Protein Name	GLYR1
Immunogen	Synthesized peptide derived from human GLYR1 AA range: 130-180
Specificity	This antibody detects endogenous levels of GLYR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Chromosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GLYR1 mouse mAb