




I17EL mouse mAb

Catalog No	YP-mAb-10991
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	IL17REL
Protein Name	I17EL
Immunogen	Synthesized peptide derived from human I17EL AA range: 27-77
Specificity	This antibody detects endogenous levels of I17EL at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	integral component of plasma membrane,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using I17EL mouse mAb