




CF218 mouse mAb

Catalog No	YP-mAb-11012
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LINC00518 C6orf218
Protein Name	CF218
Immunogen	Synthesized peptide derived from human CF218 AA range: 50-100
Specificity	This antibody detects endogenous levels of CF218 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CF218 mouse mAb