




LUR1L mouse mAb

Catalog No	YP-mAb-11020
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	LURAP1L C9orf150 HYST0841
Protein Name	LUR1L
Immunogen	Synthesized peptide derived from human LUR1L AA range: 142-192
Specificity	This antibody detects endogenous levels of LUR1L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	polymorphism:The poly-Gly region is polymorphic and the number of Gly varies in the population (from 9 to 12).,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LUR1L mouse mAb