




CG062 mouse mAb

Catalog No	YP-mAb-11023
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C7orf62
Protein Name	CG062
Immunogen	Synthesized peptide derived from human CG062 AA range: 59-109
Specificity	This antibody detects endogenous levels of CG062 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Testis-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CG062 mouse mAb