




NPPL3 mouse mAb

Catalog No	YP-mAb-11111
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NPIPL3 KIAA0220
Protein Name	NPPL3
Immunogen	Synthesized peptide derived from human NPPL3 AA range: 831-881
Specificity	This antibody detects endogenous levels of NPPL3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,similarity:Belongs to the NPIP family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NPPL3 mouse mAb