




ACPL2 mouse mAb

Catalog No	YP-mAb-11144
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ACPL2 UNQ370/PRO706
Protein Name	ACPL2
Immunogen	Synthesized peptide derived from human ACPL2 AA range: 201-251
Specificity	This antibody detects endogenous levels of ACPL2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein . Colocalizes to Golgi apparatus in a B3GAT3-dependent manner. .
Tissue Specificity	Widely expressed. Strongly expressed in spleen, fetal liver, moderately in placenta, pancreas, kidney, thymus and colon.
Function	catalytic activity:A phosphate monoester + H(2)O = an alcohol + phosphate.,similarity:Belongs to the histidine acid phosphatase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ACPL2 mouse mAb